

From: et@lmg.usap.gov
Sent: Saturday, November 05, 2016 4:47 AM
To: et@lmg.usap.gov
Subject: Oxygen System Weekly Maintenance: Sat Nov 5 04:47:01 GMT 2016

Oxygen System Weekly Maintenance

Date: Sat Nov 5 04:47:01 GMT 2016

There are clickable links in this email with detailed step by step instructions if you are unsure of any operation

TURNED SYSTEM OFF YESTERDAY
AT 68W

- N/A* 1. [Connect to the oxygen instrument using RDP](#)
- o click STOP in main window
 - o close program with (red X)
2. Copy all new data files (*.mr and *.hr from DesktopData) to the USB drive
While those are busy copying (it takes a while), continue with:
- N/A* 3. [Swap out the O2 system moisture trap](#) and note new trap positions *LEFT AS IS*
- o new trap "#", now in chiller: _____
 - o old trap "#", now in dryer: _____
4. Open the right side of the cylinder box and note the high-side and low-side regulator pressures
- o HS 700 / 3
 - o MS 1000 / 3.5
 - o LS 750 / 3
 - o LT 1800 / 3.5
 - o WT1 800 / 4.0
 - o WT2 1000 / 4.0
5. Ideally, all six low side cylinder pressures read 3.0 +/- 1.0 psig. Note any outliers
6. Close up the cylinder box
- N/A* 7. Check if the peristaltic pump inside the Manifold box is spinning. _____ (yes/no)
- N/A* 8. Confirm PC time is syncing with GMT time _____ (yes/no)
- N/A* 9. Select Run | Start in the MS VB screen
- N/A* 10. Click 'Enter in AutoRun' when prompted
- N/A* 11. Change the 'Trap in Use' to the new trap no. now in chiller
- N/A* 12. [Check that FLWTb reads 60 +/- 5 sccm, FLSP reads 100 +/- 10 sccm, and Pfridge is < 1200 torr](#)
- N/A* 13. Confirm no USB errors - note otherwise
14. Log the event in the Cruise Data Report
15. Check the fan on the 01 deck intake. Running? _____ (yes/no)
- N/A* 16. Date and Time (GMT) of weekly maint 210 5 NOV 2016
17. Scan this log sheet, store it on the D: drive of the ET computer and email a copy to Britt Stephens (stephens@ucar.edu) and Andrew Watt(watt@ucar.edu)